

easyflow®

"WE CARE YOUR HEALTH"

VIEW WEBSITE



www.ermedmedikal.com.tr







ERMED TIP MEDİKAL

We, "ERMED TIP MEDIKAL" established in 2004 as a retail business, has developed rapidly over time and become one of the leading manufacturer & wholesale supplier of medical disposables and PPE (Personal Protective Equipment) in Turkey.

With our warehouse of 9.500 m² area, located in Kocaeli/Turkey, we consistently keep ready stock to cover the demands of our customers such as more than 400 private hospitals all over Turkey and centralized public tenders mainly with our patented brand name EASYFLOW®.

With the factory of 5.000 m² area, located in Kocaeli/Kartepe,Turkey, we provide opportunity of fast delivery and competitive offers to our valuable customers by manufacturing high quality and volume of several product;



easyflow®

ORTHOPAEDICS

PLASTER OF PARIS BANDAGE

CAST PADDING

SYNTHETIC CASTING TAPE

ROLL SPLINT ELASTIC BANDAGE

SELF-ADHESIVE BANDAGE

ELASTIC NONWOVEN SURGICAL TAPE

COTTON SURGICAL TAPE SILK SURGICAL TAPE PAPER SURGICAL TAPE INJECTION PLASTER

INCONTINENCE

ADULT DIAPER WITH PANTS **ADULT DIAPER WITH TAPES**

UNDER PAD

BABY DIAPER

PAWPAD VETERINARY DIAPER

PERINEAL AREA AND BODY CLEANSING TOWEL

BARRIER SPRAY

ADHESIVE REMOVER SPRAY **BODY CLEANSING FIBER** HAIR WASHING CAP

DISPOSABLES

GAUZE BANDAGE

GAUZE DRESSING

ROLL GAUZE DRESSING

MAKEUP REMOVAL PAD

COTTON DENTAL ROLL

COTTON BUDS

COTTON PAD

NEBULIZIER MASK

OXYGEN MASK

NASAL CANNULA MANUEL RESUSCITATOR

AEROCHAMBER

VAGINAL SPECULUM

3 BALL EXERCISER ENEMA SET

SURGICAL HAND BRUSH

ALCOHOL SWAB

GEL COMPRESS

STERILIZATION ROLL

WRAP PAPER

WOODEN TONGUE DEPRESSOR

TOURNIQUET

SHARP CONTAINER DRUG BOWL

OTOSCOPE TIP

EMESIS BASIN

URAINARY

URINAL

PATIENT WRIESTBAND **ECG ELECTRODE**

PERSONAL PROTECTIVE EQUIPMENT

PROTECTIVE COVERALL PROTECTIVE COVERALL TAPERED

FOOT PROTECTOR

FACE SHIELD

SURGICAL MASK

SURGICAL HAIR BONNET SHOE COVER

ANTI EMBOLISM STOCKING

PATIENT GOWN


SURGICAL GOWN

PAGE 04-07

PAGE 08-13

PAGE 14-19

PAGE 20-23





PLASTER OF PARIS BANDAGE

Shaped easily due to its soft composition.

Areas of usage; plastering applied by circular rotations without fixing in fractures, orthopedic corrections, joint and bone irregularities.

Made of 100% cotton yarn – gauze bandage and orthopedic plaster, and the time for starting for hardening takes place between; 2-8 minutes.

Cast with quick setting time model is available.

★ Oem available / Different sizes available.

REF NO: ESY3.01000

SIZES			
5cm 7,5cm 10cm 15cm 20cm			



SYNTHENTIC CASTING TAPE

EASYFLOW Is a synthetic fill-band which is impregnated with polyurethane resin, which is activated by water prior to use.

Excellent Molding - provided excellent multi- faceted stretching and excellent molding in curved area.

Smooth Coverage - Top quality sub- layer and special resin layer provides smooth surface and clean and strong end- lamination.

Cost Effective - Compared to POP band, requires less rolls and refilling.

Fiberglas and Polyester Selection - For those who prefer non-fiberglass sub-layer, top quality polyester fill band is available.

Lightness, incredible power, endurance, air permeability and X-ray transparency.

Soft cast model made of semi-rigid polyester (non - fiberglass) which allows earlier mobility and reduces most of the problemswith rigid immobilizaon is available.

Oem available / Different sizes available.



CAST PADDING

Minimizes the formation of wounds sourced from friction and provides the patient to get healed under comfortable conditions.

It can be torn and separated easily, therefore it provides contribution for adaptation and application, and with under plaster cotton, it minimizes the exposure of weight over the fracture after padding.

It has an excellent adaptation property thus it prepares easy application and an environment for smooth plaster surface.

★ Oem available / Different sizes available.

REF NO: ESY3.02000

SIZES				
5cm 7,5cm 10cm 15cm 20cm				



SYNTHENTIC ROLL SPLINT

Fabric impregnated with very powerful 8-layered polyurethane resin.

Fill which does not include soft knitting in both sides.

Use whitout gloves.

Excellent molding for curved applications.

A resin basis- splint which is packed, curable with water.

Various dimensions which have been cut previously for different applica tions.

Economic roll type product with cut and close- folio packaging.

★ Oem available / Different sizes available.

REF NO: ESY1.01000

 SIZES

 2 inch - 5 cm
 3 inch - 7.5 cm
 4 inch - 10 cm
 5 inch - 12.5 cm
 6 inch - 15 cm

REF NO: ESY1.02000

SIZES				
2 inch - 5 cm	3 inch - 7,5 cm	4 inch - 10 cm	5 inch - 12,5 cm	6 inch - 15 cm

ELASTIC BANDAGE

It is made of polyester for weft yarn, whereas made of cotton yarn and latex for warp yarn. As it is made of cotton intensively, no allergic reaction occurs on the skin.

Elastic bandage may be used by everybody as a fixing instrument for used for the fixing of compression for muscular-based pains in case of slight muscular damages.

It has 125% elongation - with mainly 2 metal fixing clamp.

- $\slash\hspace{-0.4em}\cancel{x}$ Different width of bandage can be manufactured from 5cm to 20cm
- ★ Oem available / Different sizes available.













SELF- ADHESIVE BANDAGE

Cohesive bandage which is suitable for medium support or compression.

Its adhesive properties are long lasting and provide firm, controlled support.

Fluffy non-fray edges prevent the bandage from rucking or marking.

Limitation of movement (sprains & strains) and Correction of false position & posture.

Water repellent.

Air permeable.

Thinner substrate.

★ Different sizes are available.

REF NO: ESY3.03000

SIZES	
2,5 cm x 4.5 m	
5 cm x 4.5 m	
7.5 cm x 4.5 m	
10 cm x 4.5 m	
15 cm x 4.5 m	



ELASTIC NONWOVEN SURGICAL TAPE

Nonwoven surgical tape with hypoallergenic adhesive


for securing medical devices and exation of dressings.

Gentle and secure adhesion for sensitive skin, hypoallergenic.

Permeable to air and moisture.

Easy, straight, bi-directional tear.

★ Different sizes are available.



REF NO: ESY3.04000



SILK SURGICAL TAPE

Surgical silk tape for securing medical

devices and exation of dressings

Silk acetate fabric carrier

Powerful adhesion, can be torn easily by hand

Strong, reliable securing

Skin friendly adhesive

★ Different sizes are available

REF NO: ESY3.06000

SIZES			
5m x 2., cm	5m x 5 cm	5m x 10 cm	10m x 5 cm

INJECTION PLASTER

Hypoallergenic

Absorbent pad that does not adhere at the wound

Breathable

The adhesive does not leave waste

Individually packed





COTTON SURGICAL TAPE

Zinc oxide surgical tape with hypoallergenic adhesive for securing medical devices and exation of dressings.

Gentle and secure adhesion for sensitive skin, hypoallergenic.

Permeable to air and moisture.

Easy, straight, bi-directional tear.

★ Different sizes are available.

REF NO: ESY3.05000

SIZES			
5m x 2,5 cm 5m x 5 cm 5m x 10 cm 10m x 5 cm			





PAPER SURGICAL TAPE

Paper surgical tape with hypoallergenic adhesive for securing medical devices and exation of dressings.

Polyester viscose non-woven carrier.

Gentle and secure adhesion for sensitive skin, hypoallergenic.

Permeable to air and moisture.

Highly conformable, also on contoured parts of the body.

♯ Different sizes are available.

REF NO: ESY3.07000

SIZES				
5m x 1,25 cm 5m x 2,5 cm 5m x 5 cm 10m x 1,25 cm				

REF NO: ESY3.08000






















ADULT DIAPER WITH TAPES

Textile Outer Surface

Special Side Strip System

Super Absorbent Particles

Extra Dry Layer

Leak-proof Barriers

Suitable for Sensitive Skin

Wetness Indicator

Waist Measurement

REF NO: ESY3.09000

SIZES			
Small	Medium	Large	X-Large



Super Absorbent Particles



Textile Backsheet



Leakproof Barriers

ADULT DIAPER WITH PANTS

Textile Outer Surface

Special Side Strip System

Super Absorbent Particles

Extra Dry Layer

Leak-proof Barriers

Suitable for Sensitive Skin

Wetness Indicator

Waist Measurement

REF NO: ESY3.10000

SIZES			
Small Medium Large X-Large			



Super Absorbent Particles



Leakproof Barriers



Leakproof Barriers



Extra Dryness Layer



Suitable For Sensitive Skins

UNDER PAD

EASYFLOW Disposable Underpad, 100% breathable and liquid Proof Backsheet. Thanks to the Super Absorbent Particles, the inner pad absorbs liquid quickly and passes it to the substrate that prevents the mattress from getting wet and getting dirty. The sealed outer surface, which does not let the liquid out, prevents slipping. Thanks to the availability in various sizes, the underlay can be used in a variety of ways, for example to protect beds / mattresses, chairs / seats and other furniture and surfaces.

Cotton-like Top Layer

Super Absorbent Texture

Leak-proof Bottom Layer

Special Slip Preventing System

No skin damage

10-30 pieces in 1 pack

REF NO: ESY1.04000

	SIZES	
60 x 60 cm	60 x 90 cm	90 x 180 cm









DISPOSABLE CHANGING MATS

Disposable changing mats are highly convenient for use both at home and during travel.

Super Easy to change Baby Diaper, no hassle of washing and cleaning.

Non-woven absorbent material

Suitable for newborns and their sensitive skin

Easy to use and disposable after usage

Easy to use and disposable after usage

★ Hygienic Diaper Changing (can also be used as bed/chair protection)



REF NO: ESY1.05000

	SIZES	
60 x 60 cm	60 x 90 cm	90 x 180 cm







Comfort Friendly



Soft Tissue

PUPPY TRAINING PAD



It is used for animals.

It has good absorption.

Inner pad absorbs liguid quickly.

The sealed outer surface, which does not let the liquid out, prevents slipping.



Protective Layer

REF NO: ESY1.06000

	SIZES	
60 x 60 cm	60 x 90 cm	90 x 180 cm



Paw-Friendly

PERINEAL AREA AND BODY CLEANSING TOWEL

REF NO: ESY3.12000













Easyflow Perineal and Body Cleansing Towel are designed to clean the whole body in cases where is not posible to take a shower. Packaged ready to use, Easyflow Perineal and Body Cleansing Towel cleanse, Moisturize, nourish and protect the skin. Compatible with skin pH. Does not contain alcohol and does not require extra rinsing or drying. It helps you to have a healthier skin with its rich content. It has the feature of preventing and removing bad ador ind the perineal and other areas. It leaves a fresh and clean feeling after use.

Usage Instruction

Open the package from the part shown on the label. Apply to the area you want to apply by wiping. Use a separate towel for each zone.

Warnings and First Aid Precautions:

Do not use if package is damaged. Do not use on open wounds. Do not keep the product on the wound, do not use it as a buffer. Do not reuse. For 0-6 months babies, breastfeeding or pregnant women should be used under the supervision of a doctor. Do not use other cleaning products before, after or with this product. Avoid contact with eyes. If contact occurs, wash thoroughly with water. For external use. Keep out of reach of children. Keep away from sunlight. Store at room temperature.

REF NO: ESY3.11000





BARRIER SPRAY

Silicon - based item.

Provides skin - moisture.

Accelerates the treatment process.

Creates barrier against body liquids.

Obstructs liquid contact.

Packed in single - boxes.

REF NO: ESY3.28000



ADHESIVE REMOVER SPRAY

Silicon - based item.

Provides plaster etc. adhesives removed easily.

No skin damage.

Cohorent with sensitive skins.

Sold in single - boxes.

Adhesive remover hankies may be given along with.

REF NO: ESY3.13000

BODY CLEANSING FIBER

Body Cleansing Bag, offers a highly effective and easy care opportunity thanks to its hand-wearablebag-shaped structure. With its honeycomb texture, it easily removes even bulky remains from the body. With its 85 GSM thick soft cloth material, it contains more solution than wipes. One package of Body Wash Bag canbe used for a person's entire body Care can be provided. It does not dry the skin as in cleaning processes using soapand water. It gives the skin elasticity and contributes to healthier skin. It allows body cleansing to be done inless time compared to traditional washing methods.

Usage:

Open the package from the marked place. Add enough water on the fiber to cause foaming After cleaning, dry the area with a dry and clean cloth.

Cautions:

Keep away from children, do not expose to direct sunlight.

Store in dry and cool place (at room temperature)

It should not be used on open wounds, genital area and allergic skin

For external use



HAIR WASHING CAP

It and provides cannot be hair and washed. scalp It is care a and disposable cleaning hair for washing people, who bonnet are that does bedridden, not have require limited, rinsing mobility, and contains have special shampoo needs conditione.

Use with Microwave;

Keeping the opening label of the package half open, place it in the microwave oven. Soak for 20 seconds at max 700 W.Check the temperature of the product before use. Do not use at temperature above 35-40 C Open.





Cilt Dostu Skin Friendly





REF NO: ESY3.14000







GAUZE COMPRESS

Made of hydrophilic gauze bandage produced from 100% cotton yarn, having higher permeability, having no dye.

Composed of tight and regular knitted gauze bandage.

No yarn dangles from any part of the sponges.

4-8-12-16 layer options are available and may be used in various areas.

X-Ray (Radiopaque) option is available. Radiopague X-Ray is blue in color, and is resistant against vapor autoclave sterilization. As the sponge can be easily determined to be having radio-pact property (XRY) in the X-Ray device, it can be used in surgical operations.

Gauze compress used in abdomen is available

SIZES (FLAT-SURGICAL)

4 LAYERS

Customized size of compress can be manufactured

5 x 5 cm

Compress can be manufactured as non-sterile or sterile

REF NO: ESY1.07000

7,5 x 7,5 cm

SIZES (FLAT-SURGICAL)

4 LAYERS

8 LAYERS	5 x 5 cm	8 LAYERS	7,5 x 7,5 cm	
10 LAYERS	10 LAYERS 5 x 5 cm		7,5 x 7,5 cm	
12 LAYERS	5 x 5 cm	12 LAYERS	7,5 x 7,5 cm	
14 LAYERS	5 x 5 cm	14 LAYERS	7,5 x 7,5 cm	
SIZES (EL AT	T-SURGICAL)	SIZES (FLAT-SURGICAL)		
OIZEO (I EA	i dontoidat,	·	,	
4 LAYERS	10 x 10 cm	4 LAYERS	8 x 18 cm	
·		•		
4 LAYERS	10 x 10 cm	4 LAYERS	8 x 18 cm	
4 LAYERS 8 LAYERS	10 x 10 cm	4 LAYERS 8 LAYERS	8 x 18 cm 8 x 18 cm	



GAUZE DRESSING

It is made of Viscose/Cotton for weft yarn, whereas made of polyamide yarn for warp yarn.

As its both sides are woven, no yarn is left from the both sides, and it can be used as sterile.

Elastic ratio is; 100%

Special sizes can be produced.

Elastic type gauze dressing is available, which is to be used for the cases in which slight compression is necessary.

Different with of bandage can be manufactured from 5cm to 20cm

Different length of bandage can be manufactured



ROLL GAUZE DRESS

It is the roll bandage made of 20 wire cotton yarn.

It is mainly 100 meter length roll gauze dressing.

Various sizes are available.

Length: 10-100 meters

Different with of dressing can be manufactured.

X-ray (Radiopaque) line dressing is available.



COTTON

Hydrophilic roll cotton is used in hospital and polyclinics as a personal hygiene or daily use.

Hydrophilic roll cotton is a fast absorber which enables water retention.

Absorption occurs within 5 seconds.

50 gr 100 gr 500 gr 1000 gr

REF NO: ESY3.18000



REF NO:ESY26000

COTTON SWABS

%100 natural cotton

Beauty Care Apply eyeshadow, remove makeup and nail polish.

Baby Care Reach and softly clean a baby's delicate skin.

First Aid Perfect for general first aid to apply aintments and creams.

Cleaning Perfect for Reach those small, hard-to-clean areas around the house.





REF NO: ESY3.16000

MAKEUP REMOVER PAD

100% cotton

Soft texture and does not lint.

Perfect for applying lotion and removing make-up.

REF NO: ESY1.10000





COTTON PAD

Cottonpad; is formed in such a manner that; it has been regularly expanded inside the cotton hydrophilic gauze bandage, folded entirely by gauze hydrophilic gauze, whereas no gibers and cotton hang out.

The cotton which is required to take place within the cotton- pad is minimum $\,10\,\mathrm{gr.}$

Woven by 100% cotton yarn and made of tight and woven gauze bandage. $\,$

10cm x 15cm, 10cm x 20cm, 15cm x 20cm ve 12cm x 24cm sizes variety.



REF NO: ESY3.19000

REF NO: ESY3.23000



NEBULIZER MASK

Made of medical grade pvc.

Non-sterile.

8ml capacity - drug chamber.

Ergonomic mask appropriate with features.

Channeled and resistant against breakage and bending 200cm long connection hose.

Connector consistent with standards.

Disposable.

AEROCHAMBER

Lower pressure design provides continuous and gentle vacuum and minimizes the tissue trauma.

Volume marks enable easy measuring.

To place it easily over the surfaces, bottom part has been designed in plane shape.

Hanger for safe binding and bed side clip.

REF NO: ESY3.20000







OXYGEN MASK

Made of medical grade pvc.

Non-sterile.

8ml capacity – drug chamber.

Ergonomic mask appropriate with features.

Channeled and resistant against breakage and bending 200cm long connection hose.

Connector consistent with standards.

Disposable.

VAGINAL SPECULUM

It is used for vaginal examination, single use.

Smooth and slippery surface.

No sharp edges.

Screwed and standard models are available.

SIZES			
Small	Medium	Large	



RFF NO: FSY3 21000



REF NO: ESY3.25000

NASAL CANULLA

Channeled and resistant against breakage and bending 200cm long connection hose.

Connector consistent with the standards.

Disposable.

3 BALLS EXERCISER

Using a respiratory exerciser after surgery will help you keep your lungs clear.

The device will also help keep your lungs active when you are recovering from surgery as if you were at home performing your daily activities.

It will help you breathe properly to inhale and exhale morea.

It is used to reduce breathing problems and shortness of breath.

REF NO: ESY3.22000



REF NO: ESY3.26000



MANUEL RESUSCITATOR

It is a manually applied breathing apparatus to the person who wants to have artificial respiration.

PVC and silicone options are available.

Adult, pediatric and neonatal models are available.

The package includes an oxygen tube, reservoir, mask, berman airway, mouth opener.

SIZES				
Adult	Pediatric	Neonatal		

ENEMA SET

Rectum end available for easy use.

1,750 ml capacity.

Made of soft PVC material.

Hose resistant to 110cm bending Castile soap.


Liquid tight napkin.

PE glove.

Offset gel.

SURGICAL HAND BRUSH

REF NO: ESY3.27000





Impregnated with Chlorhexidine Gloconate (4%) or Povidine-Iodine. (7.5%)

It has ergonomic design.

Soft bristles for no skin irritation.

Smooth and soft sponge absorbs up to 90 ml liquid.

REF NO: ESY3.31000

STERILIZATION ROLL



Non-toxic, highly resistant, reinforced film

Medical grade kraft paper with a sterilization barrier, available in 60 gsm and 70 gsm.

Can be easily sealed with all types of sealing devices.

Clear and precise indicator color changes.

Opens without tearing or leaving fibers.

★TYVEK STERILIZATION ROLL available.

REF NO: ESY3.32000

WRAP PAPER



Consistent with all international standards.

Permeable and forms bacterial barrier.

Produced from pure cellulose.

Colors: Green, Blue, White.

REF NO: ESY3.29000

REF NO: ESY3.33000





ALCOHOL SWAB

Alcohol swabs are used by people who use drugs to clean an injection site before injection.

It is also used as to clean their fingers and thumb before an injection and to remove any blood from the injection on their fingers and other surfaces.

It helps ensure reliable site preparation for safe, hygenic insulin injections.

Includes 70% isopropyl alcohol.

WOODEN TONGUE DEPRESSOR

Made of Class 1 white birch.

Soft and smooth surface.

Stored in homogeneous -100 packages.

Sterile and non-sterile versions.





GEL COMPRESS

They are products that are used cold when left in the freezer, and used as hot when heated and serve as heating.

It does not freeze in the freezer and maintains its coldness for a long time.

It maintains its temperature for a long time when placed $% \left(\mathbf{r}\right) =\left(\mathbf{r}\right)$

in a microwave oven or in hot water.

It is used with its plastic packaging.

Offered in 3 different sizes.

TOURNIQUET

Applied on bones. (Arm-biceps, leg-femur)

Tourniquet should be applied for serious bleeding cases which are unable to be prevented by other procedures and having a serious danger of life.

Tourniquet is necessary in cases in which large veins (artery) are cut.

Tourniquet should be placed slightly over the wound, not to match with joint and without touching the wound.

REF NO: ESY3.35000 REF NO: ESY3.39000





SHARP CONTAINER

Used for medical waste.

Quality plastic material.

Lockable cover.

SIZES					
1,2lt	2lt	3lt	5lt	10lt	30lt

URINARY

Used for patients' need to go to bathroom.

Patient urinary has cover and is made of plastic.

Washable.

Used for the patients who are unable to leave their beds especially.

Appropriate for use by male and female.

REF NO: ESY3.36000







DRUG BOWL

Made of firm plastic.

Includes cover for drug bowl.


URINAL

High quality plastic urinal (bed shovel) made of transparent plastic.

Light and easy to use.

Washable.

REF NO: ESY3.37000



OTOSCOPE TIP

Disposable

Produced from materials which can be re-sterilized.

1 no: 2 mm 2 no: 3 mm 3 no: 4 mm

REF NO: ESY3.41000

PATIENT WRISTBAND

Adult ad pediatric, White, red, green and blue.

Mother baby.

Blue and pink.

Card and no card options.

Various color preferences.

Locking system.

REF NO: ESY3.38000





EMESIS BASIN

Used for medical dressing wastes.

Disposable.

Black and white color options.

REF NO: ESY3.42000



ECG ELECTRODE

ECG electrodes are hard gel Ag-AgCI electrodes for electrocardiography.

Higher conductive gel prevents conductivity irregularities.

Provides firm adhesion on the skin even against dirt and contamination.

Due to the support of gel and skin's wider and healthy contact surface, electrical signals are transmitted without beng exposed to external factors.

Adult, paediatric, neonatal opitons are available.

















PROTECTIVE COVERALL

Safe and lightweight protection for a wide variety of uses. The protective suit meets the requirements of EN 14126: 2003 for protection against biological hazards in the highest performance

These gyms offer maximum comfort as they are designed to provide a correct fit and ease of movement.

It has an elastic head and a self-adhesive chin cap.

It is low and solid.

Cuffs and ankles are elastic.

It has a self-adhesive zipper closure.

Thumb loop (For safety, use of thumb loops is recommended when wearing double gloves.)

Colors: White

SIZES			
Small	Medium		
Large X-Large			
XX-Large	XXX-Large		



REF NO: ESY1.13000

FOOT PROTECTOR

Thanks to the serrated surface, it is non-slip and also waterproof. Light and easy to wear.

Made of soft, durable, breathable and water-resistant fabric.

Flexible wrist ties for inlet-outlet limitation and fluid protection.

Flexible calf circumference for maximum comfort and protection

To magnify protection against hazardous chemicals and dust; The base is produced in two layers.

It has with insole and withput insole options.

Colors: White

SIZES					
42	44	46	48	50	



FACE SHIELD

It prevents the transmission of viruses

It is not difficult to assembl

It is completely hygienic

It provides clear view

It is resistant to cracking and breaking



REF NO: ESY1.12000

PROTECTIVE COVERALL TAPERED

Safe and lightweight protection for a wide variety of uses. The protective suit meets the requirements of EN 14126:2003 forprotection against biological hazards in the highest performance class.

These gyms offer maximum comfort as they are designed to provide a correct fit and ease of movement.

It has an elastic head and a self-adhesive chin cap.

It is low and solid.

Cuffs and ankles are elastic.

It has a self-adhesive zipper closure.

Thumb loop (For safety, use of thumb loops is recommended when wearing double gloves.)

Colors: White

SIZES			
Small	Medium		
Large	X-Large		
XX-Large XXX-Large			



SURGICAL MASK

Disposable (single use).

Hygienic.

Anti allergic.

Air permeability. (easy breathing)

Full and easy coherence with face.

Elastic and have lacings, 2 layered and 3 layered.

Color; Blue, White, Black, Green



REF NO: ESY1.16000

SURGICAL HAIR BONNET

Disposable (single use).

Hygienic.

Anti allergic.

Air permeability. (easy breathing)

Provides hair isolation.

Color; Blue, White, Green



SHOE COVER

Disposable (single use).

It is produced to prevent contamination from feet.

Any micron thickness produced can be done.

There are 1000 pieces in 1 package.



REF NO: ESY3.43000

ANTI-EMBOLISM STOCKING

Prior to and after embolic operation, for inpatients; it is appropriate for for post-trauma treatment; it is also appropriate for the, patients who have been exposed to hydrostatic oedema, contusion and torsio. At ankle level; 18-24mm Hg, decreases towards calf and femur, improves venous system (venous pump) internal compression disorder. Improves the femoral vain blood circulation.

Under-knee and upper-knee options

SIZES					
Small	Medium	Large	X-Large	XX-Large	XXX-Large





PATIENT GOWN

Long ad short sleeved models.

Air- permeable.

No latex content.

Colors: Black, White, Green, Blue, Navy Blue

REF NO: ESY3.45000

SURGICAL GOWN

Spunbond - Meltblown - Spunbond (3 layers)

It is produced from fabric.

It canbe produced specially in desired dimensions.

It is resistant to puncture, tearing and abrasion.

It has water repellent feature.

It maintains body temperature.

The neck part provides connection.

Its arms are ribbed.

It has an adhesive part at the neck and a tie at the waist.

Water repellent property.

Water repellent feature standard size.

It offers high comfort and comfortable working feature.

Colors : Blue

Sizes : Standard



TWO PAIRS OF KNOTS



FINE STITCHING



KNITED CUFF



ABDOMEN COMPRE ER **NASAL OXYGEN CANULL** Sł XYGE ANDAGE COMPRESS REMOVER SPRAY R m GOWN **OSCOPE** ES H N H F-A **ON PLAS1** SEL Ш CARE YOUR U





• VIEW WEBSITE•



ERMED TIP MEDİKAL TİC. A.Ş.

HeadQuarter

Köseköy Mah. Yılmaz Sok. No:15 Kartepe-KOCAELİ/ TURKEY

Phone: 0 262 223 44 80 - **Fax:** 0 262 223 44 94

in fo@ermed medikal.com.tr-www.ermed medikal.com.tr

Factory

Ovacık Mah. Kullar Cad. No:185/A -KOCAELİ/ TURKEY

Phone: 0 262 223 44 80 - **Fax:** 0 262 223 44 94

info@ermedmedikal.com.tr - www.ermedmedikal.com.tr